


OsteoBiol®
by Tecnoss

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H PHYSICOCHMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE J BIOMED MATER RES APL BIOMATER, 2010 FEB; 92(2):409-19

(2) NANNMARK U, SENNERBY L THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS CLIN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70. EPUB 2008 JAN 30

(3) CRESPI R, CAPARE P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLINE E CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

(4) CASSETTA M, PEROTTI V, CALASSO S, PIATELLI A, SINJARIB, IEZZI G BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT

(5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1; 133(6):3506-13

(6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION CLIN ORAL IMPLANTS RES, 2014 SEP 15 EPUB AHEAD OF PRINT

(7) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS CLIN IMPLANT DENT RELAT RES, 2011 NOV 14, EPUB AHEAD OF PRINT

(8) CARDAROPOLI D, CARDAROPOLI G PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

(9) BARONE A, CRESPI R, NICOLI ALDINI N, FINI M, GIARDINO R, COVANI U MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS INT J ORAL MAXILLOFAC IMPLANTS, 2005, 20: 519-525

(10) BOTTINI LP, RICCI I, PIATELLI A, PEROTTI V, IEZZI G BUCCO-LINGUAL CRESTAL BONE CHANGES AROUND IMPLANTS IMMEDIATELY PLACED IN FRESH SOCKETS IN ASSOCIATION OR NOT WITH PORCINE BONE: A NON-BLINDED RANDOMIZED CONTROLLED TRIAL IN HUMANS J PERIODONTOL, 2012 OCT 29, EPUB AHEAD OF PRINT

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(2) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOL, 2005, 54: 351-362.




(3) NANNMARK U, AZARMEH I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

(4) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006, 2: 99-106

(5) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PEROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9

(6) CALVO GUIRADO JL, PARDO ZAMORA G, SAEZ YUGUERO MR
RIDGE SPLITTING TECHNIQUE IN ATROPHIC ANTERIOR MAXILLA WITH IMMEDIATE IMPLANTS, BONE REGENERATION AND IMMEDIATE TEMPORISATION: A CASE REPORT
J IR DENT ASSOC, 2007 WINTER, 53(4):187-90

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
- (1) BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- (2) SANTAGATA M, GUARINELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489
- (3) COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- (4) CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB;29(1):59-67
- (5) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICOCANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

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(2) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMM, 2011, VOL. 9, NO. 3(S), 103-107

(3) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL 2006 DEC;77(12):1984-90

(4) SCARANO A, PIATTELLI A, PEROTTI V, MANZONI L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18





(5) IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBI JA, PEROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76 EPUB 2011 NOV 2

(6) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006; 2: 99-106

(7) BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
CLIN IMPLANT DENT RELAT RES, 2014 MAY 22 EPUB AHEAD OF PRINT

(8) SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PEROTTI V, PIATTELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILOCAC SURG, 2010 AUG; 68(8):1869-73

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(1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70

(2) CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

(3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2014 MAR 1 EPUB AHEAD OF PRINT

(4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2012 JUL 12, EPUB AHEAD OF PRINT

(5) CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

(6) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATE SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRIB, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2012 JAN 26, EPUB AHEAD OF PRINT

(7) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6. EPUB 2011 DEC 12



(8) SCARANO A, PIATELLI A, PEROTTI V, MANZONI L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18

(9) CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATELLI A, PEROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E189-95

(10) COVANI U, MARCONCINI S, CRESPI R, BARONE A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
J ORAL IMPLANTOL, 2009; 35(4):189-95

(11) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

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(2) FISCHER KR, FICKL S, MARDAS N, BOZEC I, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
QUINTESSENCE INT, 2014 AUG 6 EPUB AHEAD OF PRINT

(3) FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECEDENCE TYPE DEFECTS: A CASE SERIES
QUINTESSENCE INT, 2013;44(3):243-6

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

(3) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2011 NOV 14, EPUB AHEAD OF PRINT

(4) WACHTEL H, FICKLS, HINZE M, BOLZ W, THALMAIR THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

(5) SCARANO A, PIATELLI M, CARINCI F, PERROTTI V
REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT
J OF OSSEointegration, 2009

(6) HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
ZYGMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011;1:188-197

(7) RINNA C, UNGARI C, SALTAREL A, CASSONI A, REALE G
ORBITAL FLOOR RESTORATION
J CRANIOFAC SURG, 2005, 16(6): 968-972

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICOCHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

(3) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70. EPUB 2008 JAN 30

(4) CRESPI R, CAPPARE P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLINE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

(5) FBARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN; 14(3):373-9. EPUB 2010 MAY 11

(6) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOC® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

(7) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7. EPUB 2012 JUL 12

(8) BARONE A, BORGIA V, COVANI U, RICCI M, PIATELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2014 MAR 1 EPUB AHEAD OF PRINT

(9) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKLS
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING – A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2014 SEP 15 EPUB AHEAD OF PRINT

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICOCHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6. EPUB 2011 DEC 12

(3) SCARANO A, PIATELLI A, PERROTTI V, MANZONI L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18

(4) CASSETTA M, PERROTTI V, CALASSO S, PIATELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT

(5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

(6) HSU FY, CHUEH SC, WANG YJ
MICROSPHERES OF HYDROXYAPATITE/REconstituted COLLAGEN AS SUPPORTS FOR OSTEOBLAST CELL GROWTH
BIOMATERIALS 1999, 20:1931-1936

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(1) GRIFFITH LG, NAUGHTON G
TISSUE ENGINEERING-CURRENT CHALLENGES AND EXPANDING OPPORTUNITIES
SCIENCE 2002, 295:1009-14

(2) REDDI AH
MORPHOGENESIS AND TISSUE ENGINEERING OF BONE AND CARTILAGE: INDUCTIVE SIGNALS, STEM CELLS, AND BIOMIMETIC BIOMATERIALS
TISSUE ENG 2000, 6(4):351-59

(3) NAKASHIMA N, REDDI AH
THE APPLICATION OF BONE MORPHOGENETIC PROTEINS TO DENTAL TISSUE ENGINEERING
NAT BIOTECHNOL 2003, 9:1025-32

(4) SALASZNYK RM, WILLIAMS WA, BOSKEY A, BATORSKY A, PLOPPER GE
ADHESION TO VITRONECTIN AND COLLAGEN I PROMOTES OSTEOGENIC DIFFERENTIATION OF HUMAN MESENCHYMAL STEM CELLS
J BIOMED BIOTECHNOL 2004, 1:24-34

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[1] FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H PHYSICOCHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

[2] NANNMARK U, SENNERBY L THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS CUN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70. EPUB 2008 JAN 30

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
- 1 | COVANI U, AMERI S, CRESPI R, BARONE A
PRESERVATION OF THE ALVEOLAR PROCESS WITH HUMAN ETEROLOGOUS BONE. CONSIDERATIONS ON CLINICAL OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION
ITALIAN JOURNAL OF ORAL SURGERY, 2004, VOL 3, 1: 17-23
- 2 | CASSETTA M, CALASSO S, VOZZA I, DELL'AQUILA D
REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH CYLINDRICAL SANDBLASTED AND ACID ETCHED IMPLANTS: A PILOT STUDY
EUR JOURNAL OF PROSTHODONTICS, 2005;3(1):133-144
- 3 | ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOLOGICA, 2005 JUN;54(6):351-62
- 4 | BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT JOURNAL OF MAXILLOFACIAL IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- 5 | RINNA C, UNGARI C, SALTARELLI A, CASSONI A, REALE G
ORBITAL FLOOR RESTORATION
JOURNAL OF CRANIOFACIAL SURGERY, 2005 NOV; 16(6):968-72
- 6 | BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
INT JOURNAL OF PERIODONTAL AND RESTORATIVE DENTISTRY, 2008 JUN; 28(3):283-9
- 7 | BARONE A, SANTINI S, SBORDONE L, CRESPI R, COVANI U
A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION
INT JOURNAL OF MAXILLOFACIAL IMPLANTS, 2006 JAN-FEB; 21(1):81-5
- 8 | COVANI U, BARONE A, CORNELINI R, CRESPI R
CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL
JOURNAL OF PERIODONTAL AND RESTORATIVE DENTISTRY, 2006 APR; 77(4):722-7
- 9 | ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
INT JOURNAL OF MAXILLOFACIAL IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- 10 | TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT-MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INTERNATIONAL JOURNAL OF IMMUNOPATHOLOGY AND PHARMACOLOGY, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- 11 | BARONE A, COVANI U
MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NON-VASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS
JOURNAL OF MAXILLOFACIAL SURGERY, 2007 OCT; 65(10):2039-46
- 12 | CALVO GUIRADO JL, PARDO ZAMORA G, SAEZ YUGUERO MR
RIDGE SPLITTING TECHNIQUE IN ATROPHIC ANTERIOR MAXILLA WITH IMMEDIATE IMPLANTS, BONE REGENERATION AND IMMEDIATE TEMPORISATION: A CASE REPORT
JOURNAL OF DENTAL ASSOCIATION, 2007 WINTER; 53(4):187-90
- 13 | DEL CORSO M
SOFT TISSUE RESPONSE TO PLATELET-RICH FIBRIN: CLINICAL EVIDENCES
COSMETIC DENTISTRY, 2008, 3:16-20
- 14 | BARONE A, SANTINI S, MARCONCINI S, GIACOMELLI L, GHERLINE E, COVANI U
OSTEOTOMY AND MEMBRANE ELEVATION DURING THE MAXILLARY SINUS AUGMENTATION PROCEDURE. A COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS. CONVENTIONAL ROTATIVE INSTRUMENTS
CLINICAL ORAL IMPLANTS RESEARCH, 2008 MAY; 19(5):511-5. EPUB 2008 MAR 26
- 15 | BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT JOURNAL OF PERIODONTAL AND RESTORATIVE DENTISTRY, 2008 JUN; 28(3):283-9
- 16 | BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U
XENOGRAFT VERSUS EXTRACTION ALVEOLUS FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY
JOURNAL OF PERIODONTAL AND RESTORATIVE DENTISTRY, 2008 AUG; 79(8):1370-7
- 17 | COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
JOURNAL OF PERIODONTAL AND RESTORATIVE DENTISTRY, 2006 APR; 77(4):722-7
- 18 | CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT JOURNAL OF PERIODONTAL AND RESTORATIVE DENTISTRY, 2008 OCT; 28(5):469-77
- 19 | NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH, 2008 DEC; 10(4):264-70. EPUB 2008 JAN 30
- 20 | CALVO GUIRADO JL, GOMEZ MORENO G, BARONE A, CUTANDO A, ALCARAZ BANOS M, CHIVA F, LOPEZ MARI L, GUARDIA J
MELATONIN PLUS PORCINE BONE ON DISCRETE CALCIUM DEPOSIT IMPLANT SURFACE STIMULATES OSTEOINTEGRATION IN DENTAL IMPLANTS
JOURNAL OF PINEAL RESEARCH, 2009, 47(2):164-72
- 21 | SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V
REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT
JOURNAL OF OSSEointegration, 2009; 1(1):35-40
- 22 | COVANI U, MARCONCINI S, CRESPI R, BARONE A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
JOURNAL OF ORAL IMPLANTS, 2009; 35(4):189-95
- 23 | CALVO GUIRADO JL, GOMEZ MORENO G, LOPEZ MARI L, ORTIZ RUIZ AJ, GUARDIA J
TRAUMATIC MAXILLARY SINUS ELEVATION USING THREADED BONE DILATORS FOR IMMEDIATE IMPLANTS. A THREE-YEAR CLINICAL STUDY
MEDICINA ORAL, PATOLOGIA ORAL Y CIRUGIA BUCAL, 2010 MAR 1; 15(2):E366-70
- 24 | FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICOCHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
JOURNAL OF BIOMATERIALS SCIENCE: POLYMERS AND POLYMER LETTERS, 2010 FEB; 92(2):409-19
- 25 | GRENGA PL, REALE G, COFONE C, MEDURI A, CERUTI P, GRENGA R
HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30 PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL BLOW-OUT FRACTURE
OPHTHALMIC PLASTIC AND RECONSTRUCTIVE SURGERY, 2009 MAR-APR; 25(2):123-5
- 26 | CRESPI R, CAPPARE P, GHERLINE E
DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS
JOURNAL OF PERIODONTAL AND RESTORATIVE DENTISTRY, 2009 OCT; 80(10):1616-1621
- 27 | RINNA C, REALE G, FORESTA E, MUSTAZZA MC
MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX
JOURNAL OF CRANIOFACIAL SURGERY, 2009 MAY; 20(3): 881-4

30 | CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESIONS USING A COLLAGEN MEMBRANE WITH A THE MINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

31 | NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

32 | BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PRE-HYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9

33 | CALVO GUIRADO JL, GOMEZ MORENO G, LOPEZ MARI L, GUARDIA J, MARINEZ GONZALEZ JM, TRESGUERRES IF, PAREDES SD, FUENTES BRETO L
ACTIONS OF MELATONIN MIXED WITH COLLAGENIZED PORCINE BONE VERSUS PORCINE BONE ONLY ON OSSEointegration OF DENTAL IMPLANTS
J PINEAL RES, 2010, 48:194–203

34 | SCARANO A, PIATELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

35 | ROSSI R, MORALES RS, FRASCARIA M, BENZI R, SQUADRITO N
PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A NEW IMPLANT 3D NAVIGATION SYSTEM
EUR J ESTHETIC DENT, 2010 SUMMER; 5(2):172-88

37 | SCARANO A, CARINCI F, ASSENZA B, PIATELLI M, MURMURA G, PIATELLI A
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRAFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES
CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30

38 | PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSSE IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58

39 | SANTAGATA M, GUARINIETTO L, TARTARO G
A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N°:114-9

40 | BARONE A, RICCI M, CALVO GUIRADO JL, COVANI U
BONE REMODELLING AFTER REGENERATIVE PROCEDURES AROUND IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: AN EXPERIMENTAL STUDY IN BEAGLE DOGS
CLIN ORAL IMPLANTS RES, 2011 OCT; 22(10):1131-7

41 | SCARANO A, PIATELLI A, PERROTTI V, MANZON L, IEZZI G
MAYILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR;13(1):13-18

42 | CALVO GUIRADO JL, RAMIREZ FERNANDEZ MP, NEGREI B, DELGADO RUIZ RA, MATE SANCHEZ DE VALJE, GOMEZ MORENO G
EXPERIMENTAL MODEL OF BONE RESPONSE TO COLLAGENIZED XENOGRAFTS OF PORCINE ORIGIN (OSTEOBIOL® MP3): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13

43 | RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATE SANCHEZ DE VALJE, DELGADO RUIZ RA, NEGREI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

44 | CRESPI R, CAPRA P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLINE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

45 | HINZE M, VRIELINK L, THALMAIRT, WACHTELH, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197

46 | IEZZI G, DEGIDI M, PIATELLI A, MANGANO C, SCARANO A, SHIBI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76

47 | SLOTTE C, LINDFORS N, NANNMARK U
SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE DEFECTS WITH PREHYDRATED AND COLLAGENATED PORCINE BONE AND COLLAGEN BARRIERS: CASE PRESENTATIONS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):714-23

48 | BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

53 | BARONE A, ORLANDO B, CINGANO L, MARCONCINI S, DERCHI G, COVANI U
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489

54 | CALVO GUIRADO JL, MATE SANCHEZ JE, DELGADO RUIZ L, RAMIREZ FERNANDEZ MP
CALCULATION OF BONE GRAFT VOLUME USING 3D RECONSTRUCTION SYSTEM
MED ORAL PATOL CIR BUCAL, 2011 MAR 1, 16(2):E260-4

55 | ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATELLI M, CORVINO V, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2012 SPRING;5(1):19-33

56 | FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELEGRINO G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
IMPLANT DENT, 2012 JUN;21(3):190-5

52 | CASSETTA M, RICCI I, IEZZI G, DELL'AQUILA D, PIATELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9

57 | BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

58 | CASSETTA M, RICCI I, IEZZI G, CALASSO S, PIATELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8

59 | BOTTINI LP, RICCI I, PIATELLI A, PERROTTI V, IEZZI G
BUCCO-LINGUAL CRESTAL BONE CHANGES AROUND IMPLANTS IMMEDIATELY PLACED IN FRESH SOCKETS IN ASSOCIATION OR NOT WITH PORCINE BONE: A NON-BLINDED RANDOMIZED CONTROLLED TRIAL IN HUMANS
EUR J ORAL IMPLANTOL, SUMMER;5(2):149-61

60 | BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL 9, NO. 3 (S), 103-107

61 | FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

62 | FICKLS, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

63 | SILVESTRI M, MARTEGANI P, D'AVENIA F, FARINETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR;28(2):543-9

64 | WACHTEL H, FICKLS, HINZE M, BOLZ W, THALMAIRT
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

65 | RODRIGUEZ JG, ELDBANY RM
VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION
CLIN ORAL IMPLANTS RES, 2014 MAR 1 EPUB AHEAD OF PRINT

66 | FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY
CLIN ORAL IMPLANTS RES, 2014 MAR 31 EPUB AHEAD OF PRINT

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