OsteoBiol®
by Tecnoss

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- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (3) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINIARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (4) CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42
- (7) MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR
MATERIALS, 2017 SEP 8;10(9)
- (8) CHECCHI V, FELICE P, ZUCCELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTION IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278
- (9) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13
- (10) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (11) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRA-BONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (12) ROMBOUTS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

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- (1) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPHS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (3) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVESTIG, 2017 2018 JAN;22(1):505-513
- (4) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870
- (5) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44
- (6) SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9
- (7) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (8) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (9) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (10) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (11) BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U
EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL
CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

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(1) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATELLI A. **BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP**
BIOMED RES INT, 2016;2016:4086870

(2) ESPOSITO M, BARAUSSE C, PISTILLI R, PIATELLI M, DI SIMONE S, IPPOLITO DR, FELICE P. **POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL**
INT J OF ORAL IMPLANTOL, 2019;12(1):39-54

(3) FELICE P, BARAUSSE C, BARONE A, ZUCCHELLI G, PIATELLI M, PISTILLI R, IPPOLITO DR, SIMION M. **INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP**
INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

(4) BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U. **EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY**
J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447

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- (1) DI CARLO R, ZARA S, VENTRELLA A, SIANI G, DA ROS T, IEZZI G, CATALDI A, FONTANA A
COVALENT DECORATION OF CORTICAL MEMBRANES WITH GRAPHENE OXIDE AS A SUBSTRATE FOR DENTAL PULP STEM CELLS
NANOMATERIALS, 2019;9:604
- (2) CABALLÉ-SERRANO J, MUNAR-FRAU A, DELGADO L, PÉREZ R, HERNÁNDEZ-ALFARO F
PHYSICO-CHEMICAL CHARACTERIZATION OF BARRIER MEMBRANES FOR BONE REGENERATION
J MECH BEHAV BIOMED, 2019;97:13-20
- (3) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (4) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58
- (5) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAPH COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13
- (6) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (7) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85
- (8) POLIS-YANES C, CADENAS-SEBASTIÁN C, GUAL-VAQUÉS P, AYUSO-MONTERO R, MARÍ-ROIG A, LÓPEZ-LÓPEZ J
GUIDED BONE REGENERATION OF AN ATROPHIC MAXILLA USING HETEROLOGOUS CORTICAL LAMINA
CASE REP DENT, 2019; 5216362
- (9) HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011;1:188-197
- (10) SCARANO A, MURMURA G, MASTRANGELO F, LORUSSO F, GRECO LUCCHINA A, CARINCI F
A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592
- (11) ROSSI R, FOCE E, SCOLAVINO S
THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT
J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41
- (12) RINNA C, REALE G, FORESTA E, MUSTAZZA MC
MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX
J CRANIOFAC SURG, 2009 MAY; 20(3):881-4

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- (1) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOL, 2005 JUN;54(6):351-62
- (2) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006;2: 99-106
- (3) NANNMARK U, AZARMMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (4) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (5) SANTAGATA M, GUARINIELLO L, TARTARO G
A MODIFIED EDENTULOUS EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N.:114-9
- (6) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

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- (1) BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- (2) SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010; 36(6):485-489
- (3) COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- (4) CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67
- (5) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (6) LORENZON G, BUTTARELLO GM, CHESSA G
CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION
DENTISTRY, 2015, 5:2
- (7) DEVELIOĞLU H, OZCAN G, GULTEKIN SE, SENGUVEN B, YILDIRIM A
THE SHORT-TERM EFFECTS OF VARIOUS XENOGRRAFTS ON BONE HEALING IN RATS CRANIAL DEFECTS
BIOMEDICAL RESEARCH, 2018;29(8); 1598-1602

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- (1) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- (2) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC; 77(12):1984-90
- (3) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMM, 2011, VOL. 9, NO. 3(S), 103-107
- (4) CAKIR M, KARACA IR, AYSEGÜL F, KAYMAZ F, BOZKAYA S
EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION
CLIN ORAL IMPLANTS RES, 2015 MAR-APR; 30(2):279-85
- (5) KOLMAS J, SZWAJA M, KOŁODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENIC BIOMATERIALS USED AS BONE SUBSTITUTES
J PHARM BIOMED ANAL, 2012 MAR 5; 61:136-41
- (6) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U
CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAFTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL
J CLIN PERIODONTOL, 2017 FEB; 44(2):204-214
- (7) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006; 2: 99-106
- (8) BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
CLIN IMPLANT DENT RELAT RES, 2015 DEC; 17(6):1114-26
- (9) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U
VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTER RANDOMIZED CLINICAL TRIAL
CLIN IMPLANT DENT RELAT RES, 2015 DEC; 17(6):1114-26
- (10) IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC; 23(12):1369-76
- (11) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRIC EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (12) MARCONCINI S, GIAMMARINARO E, DERCHI G, ALFONSI F, COVANI U, BARONE A
CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL
CLIN IMPL DENT RELAT RES, 2018 Dec; 20(6):906-914

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- (1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (2) KILINC A, ATAOL M
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY
BMC ORAL HEALTH, 2017 OCT 5;17(1):126
- (3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (5) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICRO TOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVEST, 2017 2018 JAN;22(1):505-513
- (6) ESPOSITO M, BARAUSSE C, BONIFAZI L, PIATTELLI M, PISTILLI R, FERRI A, FELICE P
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VERSUS SOCKET PRESERVATION AND DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN THE MOLAR REGION: 5-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
CLINICAL TRIALS IN DENTISTRY, 2020;02(4):21-37
- (7) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6
- (8) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (9) CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- (10) BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U
CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY
CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112
- (11) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (12) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

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(4) TALLARICO M, XHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM
SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(3):263-275

(4) TALLARICO M, XHANARI E, PISANO M, GATTI F, MELONI SM
MOLAR REPLACEMENT WITH 7 MM-WIDE DIAMETER IMPLANTS: TO PLACE THE IMPLANT IMMEDIATELY OR TO WAIT 4 MONTHS AFTER SOCKET PRESERVATION? 1 YEAR AFTER LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUROPEAN JOURNAL OF ORAL IMPLANTOLOGY, 2017;10(2):169-178

(5) FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
QUINTESSENCE INT, 2014 NOV-DEC;45(10):853-60

(6) FISCHER K R, TESTORI T, WACHTEL H, MÜHELMANN S, HAPPE A, DEL FABBRO MDEL FABBRO M
SOFT TISSUE AUGMENTATION APPLYING A COLLAGENATED PORCINE DERMAL MATRIX DURING SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER CASE SERIES
CLIN IMPLANT DENT RELAT RES, 2019;1-8

(7) FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES
QUINTESSENCE INT, 2013;44(3):243-6

(8) MATOH U, PETELIN M, GASPERIC R
SPLIT-MOUTH COMPARISON OF CORONALLY ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR COLLAGEN MATRIX FOR TREATMENT OF ISOLATED GINGIVAL RECESSIONS
INT J PERIODONTICS RESTORATIVE DENT, 2019;39(3):439-446

(9) VERARDI S, ORSINI M, LOMBARDI T, AUSENDA F, TESTORI T, PULICI A, OREGIA F, VALENTE NA, STACCHI C
COMPARISON BETWEEN TWO DIFFERENT TECHNIQUES FOR PERI-IMPLANT SOFT TISSUE AUGMENTATION: PORCINE DERMAL MATRIX GRAFT VS. TENTING SCREW
J PERIODONTOL; 2019 DEC 19. EPUB AHEAD OF PRINT

ARCURI C, et al. **CLINICAL AND HISTOLOGIC STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS**. MINERVA STOMATOL, 2005 JUN;54(6):351-62

BARONE A, et al. **MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS**. INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25

RINNA C, et al. **ORBITAL FLOOR RESTORATION**. J CRANIOFAC SURG, 2005 NOV; 16(6):968-72

BARONE A, et al. **IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS**. EUR J IMPLANT PROSTHODONTICS, 2006;2: 99-106

BARONE A, et al. **A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION**. INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5

COVANI U, et al. **CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL**. J PERIODONTOL, 2006 APR;77(4):722-7

ORSINI G, et al. **HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**. J PERIODONTOL, 2006 DEC;77(12):1984-90

BARONE A, et al. **MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NON- VASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS**. J ORAL MAXILLOFAC SURG, 2007 OCT;65(10):2039-46

DEL CORSO M. **SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN: CLINICAL EVIDENCES** COSMETIC DENTISTRY, 2008, 3:16-20

BARONE A, et al. **IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES**. INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9

BARONE A, et al. **XENOGRFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY**. J PERIODONTOL, 2008 AUG;79(8):1370-7

COVANI U, et al. **BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH**. INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6

CARDAROPOLI D, et al. **PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY**. INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

CRESPI R, et al. **DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS**. J PERIODONTOL, 2009 OCT; 80(10):1616-1621

RINNA C, et al. **MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX**. J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4

CARDAROPOLI D, et al. **HEALING OF GINGIVAL RESSIONS USING A COLLAGEN MEMBRANE WITH A THE MINERALIZED XENOGRFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL**. INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

SCARANO A, et al. **PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION**. J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

SCARANO A, et al. **VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRFT WITHOUT MINISCREWS AND MINIPATES: CASE SERIES**. CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30

PAGLIANI L, et al. **A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY**. CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58

SCARANO A, et al. **MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRIC EVALUATION AFTER 4 AND 6 MONTHS**. CLIN IMPLANT DENT RELAT RES, 2011 MAR;13(1):13-18

CRESPI R, et al. **CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO**. INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

HINZE M, et al. **ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY**. ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197

IEZZI G, et al. **COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS**. CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76

SLOTTE C, et al. **SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE DEFECTS WITH PREHYDRATED AND COLLAGENATED PORCINE BONE AND COLLAGEN BARRIERS: CASE PRESENTATIONS**. CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):714-23

BARONE A, et al. **A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL**. CLIN ORAL IMPLANTS RES, 2013 JAN; 24(1):1-6

SANTAGATA M, et al. **IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS**. J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489

FESTA VM, et al. **PORCINE-DERIVED XENOGRFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS**. CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13

RAMIREZ FERNANDEZ MP, et al. **ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRFTS USED IN MAXILLARY SINUS FLOOR ELEVATION**. CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

CASSETTA M, et al. **RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**. INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9

BARONE A, et al. **A RANDOMIZED CLINICAL TRIAL TO EVALUATE AND COMPARE IMPLANTS PLACED IN AUGMENTED VS. NON-AUGMENTED EXTRACTION SOCKETS. A 3-YEAR EVALUATION**. J PERIODONTOL, 2012 JUL;83(7):836-46

BARONE A, et al. **TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING**. CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

CASSETTA M, et al. **USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES**. INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8

FICKL S, et al. **PORCINE DERMAL MATRIX FOR COVERING OF RESECTION TYPE DEFECTS: A CASE SERIES**. QUINTESSENCE INT, 2013;44(3):243-6

SILVESTRI M, et al. **SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL**. INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9

WACHTEL H, et al. **THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES**. INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

RODRIGUEZ JG, et al. **VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION**. INT J ORAL MAXILLOFAC SURG, 2013 SEP;42(9):1060-6

TRAINI T, et al. **REGENERATION OF HUMAN BONE USING DIFFERENT BONE SUBSTITUTE BIOMATERIALS**. CLIN IMPLANT DENT RELAT RES, 2015 FEB;17(1):150-62

BARONE A, et al. **FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL**. CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

BARONE A, et al. **BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY**. CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):823-30

BARONE A, et al. **THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT**. CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26

CASSETTA M, et al. **BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**. CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4

GHENO E, et al. **THE EFFECTIVENESS OF THE USE OF XENOGENEIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES**. J OSSEointegration 2014;6(2):37-42

FALISI G, GALLI M, et al. **USE OF 3D CARTILAGE SCAFFOLDS FOR THE STABILIZATION OF IMPLANTS AND BONE REGENERATION WITH THE FIT-LOCK TECHNIQUE**. ACTA ODONTOL LATINOAM 2013;26(3):167-172

FISCHER KR, et al. **STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS**. QUINTESSENCE INT, 2014 NOV-DEC; 45(10):853-60

SCARANO A, et al. **EXPANSION OF THE ALVEOLAR BONE CREST WITH ULTRASONIC SURGERY DEVICE: CLINICAL STUDY IN MANDIBLE**. INT J IMMUNOPATHOL PHARMACOL, 2011 APR-JUN; 24(2 SUPPL):71-5

SCARANO A, et al. **DELAYED EXPANSION OF THE ATROPHIC MANDIBLE BY ULTRASONIC SURGERY: A CLINICAL AND HISTOLOGIC CASE SERIES**. INT J ORAL MAXILLOFAC IMPLANTS, 2015 JAN-FEB;30(1):144-9

LORENZON G, et al. **CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION**. DENTISTRY, 2015, 5:2

THALMAIR T, et al. **DIMENSIONAL ALTERATIONS OF EXTRACTION SITES AFTER DIFFERENT ALVEOLAR RIDGE PRESERVATION TECHNIQUES – A VOLUMETRIC STUDY**. J CLIN PERIODONTOL, 2013 JUL;40(7):721-7

BARONE A, et al. **VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL**. CLIN ORAL IMPLANTS RES, 2016 NOV;27(11):E105-E115

ESPOSITO M, et al. **THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL**. EUR J ORAL IMPLANTOL, 2015;8(3):233-244

OZEL B, et al. **A NEW OPTION FOR THE RECONSTRUCTION OF ORBITAL FLOOR DEFECTS WITH HETEROLOGOUS CORTICAL BONE**. J CRANIOIMAXILLOFAC SURG, 2015 OCT;43(8):1583-8

CORBELLA S, et al. **HISTOMORPHOMETRIC OUTCOMES AFTER LATERAL SINUS FLOOR ELEVATION PROCEDURE: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS**. CLIN ORAL IMPLANTS RES, 2016 SEP;27(9):1106-22

BARONE A, et al. **CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY**. CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112.

ROSSI R, et al. **THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY**. MINERVA STOMATOL, 2016 OCT;65(5):257-68

BARONE A, et al. **ESTHETIC OUTCOME OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS BY CLINICIANS WITH OR WITHOUT EXPERIENCE: A MEDIUM-TERM RETROSPECTIVE EVALUATION**. INT J ORAL MAXILLOFAC IMPLANTS, 2016;31(6)

TALLARICO M, et al. **SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL**. EUR J ORAL IMPLANTOL, 2016;9(3):263-275

LOPEZ MA, et al. **THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES**. J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

BARONE A, et al. **CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRFTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL**. J CLIN PERIODONTOL, 2017 FEB;44(2):204-214

ALFONSI F, et al. **MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS**. CURR PHARM BIOTECHNOL, 2017;18(1):64-75

BARONE A, et al. **EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY**. J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447

SCARANO A. **TRADITIONAL POSTEXTRACTIVE IMPLANT SITE PREPARATION COMPARED WITH PRE-EXTRACTIVE INTERRADICULAR IMPLANT BED PREPARATION IN THE MANDIBULAR MOLAR REGION, USING AN ULTRASONIC DEVICE: A RANDOMIZED PILOT STUDY**. INT J ORAL MAXILLOFAC IMPLANTS, 2017 MAY/JUN;32(3):655-660

BARONE A, et al. **EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL**. CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

GIULIANI A, et al. **REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY**. CLINICAL ORAL INVESTIGATIONS, 2018 JAN;22(1):505-513

IEZZI G, et al. **MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS**. CURR PHARM BIOTECHNOL, 2017, 18, 33-44

FELICE P, et al. **INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP**. INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

CRESPI R, et al. **COMPARISON OF MAGNESIUM-ENRICHED HYDROXYAPATITE AND PORCINE BONE IN HUMAN EXTRACTION SOCKET HEALING: A HISTOLOGIC AND HISTOMORPHOMETRIC EVALUATION**. INT J ORAL MAXILLOFAC IMPLANTS, 2011 SEP-OCT;26(5):1057-62

CORBELLA S, et al. **HISTOMORPHOMETRIC RESULTS AFTER POSTEXTRACTION SOCKET HEALING WITH DIFFERENT BIOMATERIALS: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS.** INT J ORAL MAXILLOFAC IMPLANTS, 2017 SEP/OCT;32(5):1001–1017

KILINC A, et al. **HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY.** BMC ORAL HEALTH, 2017 OCT 5;17(1):126

TROIANO G, et al. **COMBINATION OF BONE GRAFT AND RESORBABLE MEMBRANE FOR ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW, META-ANALYSIS AND TRIAL SEQUENTIAL ANALYSIS.** J PERIODONTOL, 2017 SEP 12:1-17. EPUB AHEAD OF PRINT

ROSSI R, et al. **A NEW INTERPRETATION OF GUIDED IMPLANT SURGERY TO ACHIEVE AN OPTIMAL RESULT IN THE ESTHETIC ZONES.** MEDICAL RESEARCH ARCHIVES, 2017 APRIL, VOL. 5, ISSUE 4

FORABOSCO A, et al. **CONCENTRATED GROWTH FACTORS IN MAXILLARY SINUS FLOOR AUGMENTATION: A PRELIMINARY CLINICAL COMPARATIVE EVALUATION.** INT J OF GROWTH FACTORS AND STEM CELLS IN DENTISTRY, 2018;1:2-7

RAMÍREZ FERNÁNDEZ MP, et al. **COMPARISON OF TWO XENOGRIFT MATERIALS USED IN SINUS LIFT PROCEDURES: MATERIAL CHARACTERIZATION AND IN VIVO BEHAVIOR.** MATERIALS, 2017 JUN 7;10(6)

AMR AEH, et al. **XENOGENIC FLEXIBLE BONE LAMINA GRAFT: A SUCCESSFUL ALTERNATIVE TO THE AUTOGENOUS ONLY BONE BLOCK GRAFT IN ALVEOLAR RIDGE AUGMENTATION: A CLINICAL, RADIOGRAPHIC AND HISTOLOGICAL EVALUATION.** JOURNAL OF DENTAL TREATMENT AND ORAL CARE, 2017;1(1);104

ATTIA AM. **CLINICAL AND RADIOGRAPHIC EVALUATION OF CORTICO-CANCELLOUS BONE MIX XENOGRIFT (OSTEOBIOL GEN-OS) IN THE TREATMENT OF HUMAN PERIODONTAL INTRABONY DEFECTS.** EGYPTIAN DENTAL JOURNAL, 2017;63;1-10

EKSTEIN J, et al. **MARGINAL BONE LEVEL AROUND CONICAL CONNECTION TAPERED IMPLANTS WITH PLATFORM SWITCHING: A MULTICENTER RETROSPECTIVE STUDY AT 14 MONTHS FOLLOW-UP.** J OSSEOINTEGRATION, 2016 JAN-APR;8(1)

NOAMI SA, et al. **EVALUATION OF PRE-HYDRATED COLLAGENATED CORTICO-CANCELLOUS GRANULES (MP3) IN AUGMENTATION OF THE MAXILLARY SINUS (PRELIMINARY STUDY).** J OF DENT, ORAL DISORD AND THERAPY, 2014;2(3): 1-8

BARONE A, et al. **EXTRACTION SOCKET HEALING IN HUMANS AFTER RIDGE PRESERVATION TECHNIQUES: COMPARISON BETWEEN FLAPLESS AND FLAPPED PROCEDURES IN A RANDOMIZED CLINICAL TRIAL.** J PERIODONTOL, 2014 JAN;85(1):14-23

LOPEZ MA, et al. **THE USE OF VARIOUS BIOMATERIALS IN COMPUTER-GUIDED CRESTAL SINUS LIFT PROCEDURES. A REPORT ON TWO CASE STUDIES WITH VOLUME COMPARISON.** ORAL IMPLANTOL, 2016;APR-JUN 9(2)

RAMÍREZ FERNÁNDEZ MP, et al. **SEM-EDX STUDY OF THE DEGRADATION PROCESS OF TWO XENOGRIFT MATERIALS USED IN SINUS LIFT PROCEDURES.** MATERIALS, 2017;10;542

COVANI U, et al. **TISSUE STABILITY OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: A 5-YEAR PROSPECTIVE SINGLE-COHORT STUDY.** J PERIODONTOL, 2014 SEP;85(9):E323-32

KAWAKAMI S, et al. **INFLUENCE OF THE POSITION OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ASSESSED WITH CONE-BEAM COMPUTED TOMOGRAPHY: A RANDOMIZED CLINICAL TRIAL.** JOURNAL INVESTIG CLIN DENT, 2018 NOV;9(4):E12362

SCARANO A, et al. **SINUS MEMBRANE ELEVATION WITH HETEROLOGOUS CORTICAL LAMINA: A RANDOMIZED STUDY OF A NEW SURGICAL TECHNIQUE FOR MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT.** MATERIALS, 2018 AUG 17;11(8)

FIGLIUZZI MM, et al. **POSTEXTRACTION IMPLANTS IN AESTHETIC AREAS: EVALUATION OF PERIMPLANT BONE REMODELING OVER TIME.** ANN STOMATOL, 2015 MAY 18;6(1):29-34

MARCONCINI S, et al. **CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL.** CLIN IMPL DENT RELAT RES, 2018 DEC;20(6):906-914

MARCONCINI S, et al. **CLINICAL SUCCESS OF DENTAL IMPLANTS PLACED IN POSTERIOR MANDIBLE AUGMENTED WITH INTERPOSITIONAL BLOCK GRAFT: 3-YEAR RESULTS FROM A PROSPECTIVE COHORT CLINICAL STUDY.** J ORAL MAXILLOFAC SURG, 2019 FEB;77(2):289-298

BERNARDI S, et al. **SHORT VERSUS LONGER IMPLANTS IN MANDIBULAR ALVEOLAR RIDGE AUGMENTED USING OSTEOGENIC DISTRACTION: ONE-YEAR FOLLOW-UP OF A RANDOMIZED SPLIT-MOUTH TRIAL.** J ORAL IMPLANTOL, 2018 JUN;44(3):184-191

SCARANO A, et al. **A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE.** J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592

SENESE O, et al. **RETROSPECTIVE ANALYSIS OF 79 PATIENTS WITH ORBITAL FLOOR FRACTURE: OUTCOMES AND PATIENT-REPORTED SATISFACTION.** ARCH CRANIOFAC SURG, 2018 JUN;19(2):108-113

ESPOSITO M, et al. **POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL.** INT J OF ORAL IMPLANTOL, 2019;12(1):39-54

FELICE P, et al. **POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG X 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR POST-LOADING RESULTS FROM A WITHIN-PERSON RANDOMISED CONTROLLED TRIAL.** INT J OF ORAL IMPLANTOL, 2019;12(1):57-72

MATOH U, et al. **SPLIT-MOUTH COMPARISON OF CORONALLY ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR COLLAGEN MATRIX FOR TREATMENT OF ISOLATED GINGIVAL RECESSIONS.** INT J PERIODONTICS RESTORATIVE DENT, 2019;39(3):439-446

KAWAKAMI S, LANG NP, et al. **INFLUENCE OF THE HEIGHT OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ASSESSED BY CONE BEAM COMPUTED TOMOGRAPHY- A RANDOMIZED CLINICAL TRIAL.** INT J ORAL MAXILLOFAC IMPLANTS, 2019;34(1):223-232

FISCHER K R, et al. **SOFT TISSUE AUGMENTATION APPLYING A COLLAGENATED PORCINE DERMAL MATRIX DURING SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER CASE SERIES.** CLIN IMPLANT DENT RELAT RES, 2019;1-8

HIROTA A, et al. **TOMOGRAPHIC EVALUATION OF THE INFLUENCE OF THE PLACEMENT OF A COLLAGEN MEMBRANE SUBJACENT TO THE SINUS MUCOSA DURING MAXILLARY SINUS FLOOR AUGMENTATION: A RANDOMIZED CLINICAL TRIAL.** INT J IMPLANT DENT, 2019; 5(1):31

TANAKA K, et al. **SINUS FLOOR ELEVATION AND ANTROSTOMY HEALING: A HISTOMORPHOMETRIC CLINICAL STUDY IN HUMANS.** IMPLANT DENT, 2019; AUG 14

RAMANAUSKAITE A, et al. **DENTAL IMPLANT OUTCOMES IN GRAFTED SOCKETS: A SYSTEMATIC REVIEW AND META-ANALYSIS.** J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E8

FARIA-ALMEIDA R, et al. **EXTRACTION SOCKET PRESERVATION WITH OR WITHOUT MEMBRANES, SOFT TISSUE INFLUENCE ON POST EXTRACTION ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW.** J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E5

ADILGULU S, et al. **CLINICAL AND HISTOPATHOLOGICAL COMPARATIVE STUDY OF TWO EQUINE-DERIVED BONE GRAFT: A HUMAN STUDY.** J PAK MED ASSOC., 2019; 69: 1617

ROSSI R, et al. **RECONSTRUCTION OF A HORIZONTAL AND VERTICAL BONE DEFECT USING THE CORTICAL LAMINA TECHNIQUE.** MED. RES. ARCH, 2019; VOL 7 (11)

ROSSI R, et al. **HORIZONTAL RIDGE AUGMENTATION WITH THE CORTICAL LAMINA TECHNIQUE: A CASE REPORT.** MODERN RESEARCH IN DENTISTRY, 2019; 4(4)

FELICE P, et al. **IMMEDIATE, EARLY (6 WEEKS) AND DELAYED (4 MONTHS) SINGLE POST-EXTRACTION IMPLANTS: 3-YEAR POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL.** CLINICAL TRIALS IN DENTISTRY, 2019;01(1):5-23

BARAUSSE C, et al. **POSTERIOR JAW REHABILITATION USING PARTIAL PROSTHESES SUPPORTED BY IMPLANTS 4.0 X 4.0 MM OR LONGER: THREE-YEAR POST-LOADING RESULTS OF A MULTICENTRE RANDOMISED CONTROLLED TRIAL.** CLINICAL TRIALS IN DENTISTRY, 2019;01(1):25-36

VERARDI S, et al. **COMPARISON BETWEEN TWO DIFFERENT TECHNIQUES FOR PERI-IMPLANT SOFT TISSUE AUGMENTATION: PORCINE DERMAL MATRIX GRAFT VERSUS TENTING SCREW.** J PERIODONTOL, 2019 DEC 19. EPUB AHEAD OF PRINT

ASLAN S, et al. **ENTIRE PAPILLA PRESERVATION TECHNIQUE IN THE REGENERATIVE TREATMENT OF DEEP INTRABONY DEFECTS: 1-YEAR RESULTS.** J CLIN PERIODONTOL, 2017 SEP;44(9):926-932

TALLARICO et al. **MOLAR REPLACEMENT WITH 7 MM-WIDE DIAMETER IMPLANTS: TO PLACE THE IMPLANT IMMEDIATELY OR TO WAIT 4 MONTHS AFTER SOCKET PRESERVATION? 1 YEAR AFTER LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL.** EUR J ORAL IMPLANTOL, 2017;10(2):169-178

FELICE P, et al. **FOUR-MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN ATROPHIC POSTERIOR JAWS: THREE-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL.** CLINICAL TRIALS IN DENTISTRY, 2020;02(1):5-26

BALDI N, et al. **XENOGENEIC DERMAL MATRIX VERSUS AUTOLOGOUS CONNECTIVE TISSUE GRAFT VERSUS NO GRAFT AT ABUTMENT CONNECTION FOR IMPROVING AESTHETICS: 6-MONTH OUTCOMES OF A RANDOMISED CONTROLLED TRIAL.** CLINICAL TRIALS IN DENTISTRY, 2020;02(2):49-62

LUONGO R, et al. **GRAFTLESS MAXILLARY SINUS FLOOR AUGMENTATION WITH SIMULTANEOUS PORCINE BONE LAYER INSERTION: A 1- TO 5-YEAR FOLLOW-UP STUDY.** INT J ORAL MAXILLOFAC IMPLANTS, 2020 JUL/AUG;35(4):808-815

SCARANO A, et al. **CORTICAL BONE LAMINA APPROACH FOR MANDIBULAR LARGE CYSTIC DEFECT: A CASE REPORT.** J BIOL REGUL HOMEOST AGENTS, 2019 NOV-DEC;33(6 SUPPL. 2):85-91

ASLAN S, et al. **ENTIRE PAPILLA PRESERVATION TECHNIQUE: A NOVEL SURGICAL APPROACH FOR REGENERATIVE TREATMENT OF DEEP AND WIDE INTRABONY DEFECTS.** INT J PERIODONTICS RESTORATIVE DENT, 2017 MAR/APR;37(2):227-233

FELICE P, et al. **IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VERSUS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: THREE-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL.** CLINICAL TRIALS IN DENTISTRY, 2020;02(3):5-25

ESPOSITO M, et al. **IMMEDIATE LOADING OF 3 MM-DIAMETER IMPLANTS AS AN ALTERNATIVE HORIZONTAL BONE AUGMENTATION FOR PLACING 4 MM-DIAMETER IMPLANTS: ONE-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL.** CLINICAL TRIALS IN DENTISTRY, 2020;02(3):61-76

ROSSI R, et al. **CORTICAL LAMINA: A NEW DEVICE FOR THE TREATMENT OF MODERATE AND SEVERE TRIDIMENSIONAL BONE AND SOFT TISSUE DEFECTS.** INT J ESTHET DENT, 2020;15(4):454-473

IEZZI G, et al. **SPLIT CREST TECHNIQUE FOR IMPLANT TREATMENT OF AGENESIS OF THE UPPER LATERAL INCISORS: RESULTS OF A RANDOMIZED PILOT HISTOLOGICAL AND CLINICAL STUDY AT 24-MONTH FOLLOW-UP.** BRAZ ORAL RES, 2020 SEP 4;34:E118

ESPOSITO M, et al. **WIDE DIAMETER IMMEDIATE POST-EXTRACTION IMPLANTS VERSUS SOCKET PRESERVATION AND DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN THE MOLAR REGION: 5-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL.** CLINICAL TRIALS IN DENTISTRY, 2020;02(4):21-37

TH ELASKARY A, et al. **VESTIBULAR SOCKET THERAPY: A NOVEL APPROACH FOR IMPLANT PLACEMENT IN DEFECTIVE FRESH EXTRACTION SOCKETS WITH OR WITHOUT ACTIVE SOCKET INFECTION (ONE-ARM COHORT STUDY).** INT J ORAL MAXILLOFAC IMPLANTS, 2021 JAN-FEB;36(1):146-153

IGLESIAS VELÁZQUEZ Ó, et al. **SPLIT BONE BLOCK TECHNIQUE: 4-MONTH RESULTS OF A RANDOMISED CLINICAL TRIAL COMPARING CLINICAL AND RADIOGRAPHIC OUTCOMES BETWEEN AUTOGENOUS AND XENOGENEIC CORTICAL PLATES.** INT J OF ORAL IMPLANTOLOGY, 2021;14(1):41-52

FOCE E, et al. **HORIZONTAL BONE AUGMENTATION WITH THE CORTICAL LAMINA TECHNIQUE: A CASE REPORT.** BIOMED J SCI & TECH RES 35(3)-2021.

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